

**CAST GOLD CROWN PREPARATION  
TOOTH #3**

**CERVICAL MARGIN  
MARGIN/WIDTH**

**Satisfactory**

The cervical margin width (chamfer) is ideally 1.0 mm.

**Unsatisfactory**

- a) The cervical margin width is less than 0.5 mm.
- b) The cervical margin width is greater than 1.5 mm.

**“a” and “b” can be grades**

**WALLS, TAPER AND DRAW  
AXIAL TISSUE REMOVAL**

**Satisfactory**

The axial tissue removal as measured from the height of contour is ideally 1.5 mm.

**Unsatisfactory**

- a) The axial tissue removal is less than 1.0 mm (resulting restoration would be over contoured).
- b) The axial tissue removal is greater than 2.0 mm. (excessive removal of tooth structure).

**“a” and “b” can be grades**